

Honorary Lecture Esgar 2013, Barcelona

MRI of HCC , Dr KADOYA, Japan

MRI FOR PROGRESSED HCC : TAKE HOME POINTS

T2 pour detection

- FAST Suppressed >>> HASTE
- HASTE if HCC with fat
- Combinaison of both is better

Fat deposition in HCC : diagnosis

- Chemical shift GRE imaging : ****
- Fat suppressed

GRE for dynamic contrast imaging

- 3 D fat suppressed T 1 gre sequences
- KWIC KIM Eur Radiol 2013, radial sampling